

Program Mode	Therapy Type	Pulse Rate	Output Mode/ Pulse Width	Type of Pain/ Program	Benefits	You Should Feel
TENS Output Modes						
P1	TENS	15Hz	Constant	Chronic	For temporary relief of pain associated with sore and/or aching muscles in the lower back due to strain from exercise or normal household and/or work activities. For temporary relief of pain associated with sore and/or aching muscles in the upper and lower extremities (arm and/or leg) due to strain from exercise or normal household and work activities.	Continuous comfortable tingling. The underlying pain should help decrease gradually after treatment.
P2	TENS	60Hz	Burst	Acute		Comfortable pulsing sensation. The underlying pain should decrease.
Р3	TENS	60Hz	Constant	Chronic		Comfortable pulsing sensation. The underlying pain should decrease.
P4	TENS	2-60Hz	Modulation	Chronic		Variable comfortable tingling and pulsing sensation (sensation should appear to come in waves). Pain should ease and there should be relief after treatment.
P5	TENS	60Hz	Modulation	Chronic		Variable comfortable mild tingling sensation (sensation will appear to come in waves).
P6	TENS	7-60Hz	Modulation	Chronic		Variable comfortable pulsing and pumping action (action will appear to come in waves).
P7	TENS	60Hz	Modulation	Chronic		Variable comfortable tingling and pumping action (action should appear to come in waves).
P8	TENS	2.45-245 Hz	Cycle	Arthritis	For symptomatic relief and management of chronic, intractable pain, and relief of pain associated with arthritis.	Variable comfortable tingling and pulsing sensation (sensation should appear to come in waves). Pain should ease and there should be relief after treatment.